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—

—

Distribution:  
1. Attachment to your Service Report

# System Check Out Protocol

OPMI Sensera, PROergo, Movena, Movena Dental, S7

APPL. N/A

- |                                     |                                     |                                       |   |
|-------------------------------------|-------------------------------------|---------------------------------------|---|
| <input checked="" type="checkbox"/> | <input type="checkbox"/>            | <b>Installation</b>                   | → Fill out chapter 1, 3, 4 and use either the Mobile X Customer perfect order application or the "Installation and Acceptance Report" |
| <input type="checkbox"/>            | <input checked="" type="checkbox"/> | <b>Preventive Maintenance</b>         | → Fill out chapter 2, 3 and 4   |
| <input type="checkbox"/>            | <input checked="" type="checkbox"/> | <b>Corrective Maintenance</b>         | → Fill out chapter 2, 3 and 4   |
| <input type="checkbox"/>            | <input checked="" type="checkbox"/> | <b>Safety Check</b>                   | → Fill out chapter 3 and 4  |
|                                     |                                     | or                                    | or  |
| <input type="checkbox"/>            | <input checked="" type="checkbox"/> | <b>Safety Check with local report</b> | → Fill out local Check List   |
| <input type="checkbox"/>            | <input checked="" type="checkbox"/> | <b>Further documents attached</b>     |   |

(The electrical safety test must be carried out in accordance with applicable local regulations and directives.)

<b>System Serial Number:</b>	<u>6628189108</u>	<b>Software Version:</b>	<u>na</u>
		<b>Build Number:</b>	<u>na</u>

**Note:** Any section has to be filled in by marking "NA"(not applicable) or "OK" / "Not OK". If you have to mark any check box "not OK", please explain in "Comments".

## Chapter 1

Follow SM-30-0032-A02

Product accessories used	Type	Serial number
303771-9020-000	Tubo Foldable	40000265

Basic:	N/A	OK	Not OK
OPMI Sensera	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
OPMI PROergo	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
OPMI Movena	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
OPMI Movena Dental	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Optional:	Serial Number
S7 Ceiling mount	
S7 Floor stand	6628189108
Inverter Tube	na

Function test: Check...	N/A	OK	Not OK
the optical Image over complete zoom range	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
zoom and focus movement over complete range	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
XY movement over complete range/centering	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
handgrip functions	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
if field of view is illuming homogeneously	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
retina protection	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
if system is undamaged, clean and in good condition	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

<b>Function test: Check...</b>	<b>N/A</b>	<b>OK</b>	<b>Not OK</b>
the tilt function of OPMI head	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
the DOF function	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
the eyepieces +/- setting	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
the bine tube for distance/movement	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

<b>Optional: Check...</b>	<b>N/A</b> <input type="checkbox"/>	<b>OK</b>	<b>Not OK</b>
movement left/right/tilt and fixation		<input type="checkbox"/>	<input checked="" type="checkbox"/>
the optical image in all possible steps		<input checked="" type="checkbox"/>	<input type="checkbox"/>
the focusing function		<input checked="" type="checkbox"/>	<input type="checkbox"/>
resight		<input type="checkbox"/>	<input checked="" type="checkbox"/>
mounted correctly		<input checked="" type="checkbox"/>	<input type="checkbox"/>
image		<input type="checkbox"/>	<input checked="" type="checkbox"/>

## Chapter 2

<b>Each year</b>	<b>N/A</b> <input type="checkbox"/>	<b>OK</b>	<b>Not OK</b>
Check if the system is undamaged, clean and in good condition	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Check if all manuals and identification labels are available, undamaged, complete and legible	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Check if the complete mechanical structure including all joints / bearings / moving parts / etc. is free of play and undamaged	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Check if the additional parts and accessories are attached properly and without play	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Check the function of the XY coupling (movement over full range/ XYZ reset / calibration...)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Check the motorized and manual tilt function of the OPMI head	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Check the Resight function (vision / focus / in – out / selftest)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Check the Invertertube function (flip)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Check if mounted accessories / eyecups and Asepsis caps are complete and in good condition	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

<b>Every 2<sup>nd</sup> year</b>	<b>N/A</b> <input type="checkbox"/>	<b>OK</b>	<b>Not OK</b>
Check if there are any signs of crack formation or wear on the supporting structure	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Check if there are any noticeable cracks, signs of wear or damage on connecting elements, suspensions and couplings and the mount of the compensating weights	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Check the camera head for function / adjustment	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

<b>Product accessories used</b>	<b>Type</b>	<b>Serial number</b>

## Chapter 3

- Reference "Installation and Acceptance Report"
- In addition follow your national regulations and if applicable all national standards!

	CZ Service	Customer*
The Safety test according to EN/IEC 62353 is performed by:	<input checked="" type="checkbox"/>	<input type="checkbox"/>

\* In case the customer is responsible, no further action is required from FSE, but it is mandatory to add this information to your service report.

Visual inspection (EN/IEC 62353)	N/A	OK	Not OK
-------------------------------------	-----	----	--------

- All marking, labels and labeling are present, complete and legible
- All protective earth connections and power cords are undamaged
- All externally accessible fuses comply with the manufacturer's specification
- The mechanical parts are undamaged
- No damage or contamination that may impair safety is present
- The associated accessories are tested together with the ME equipment /ME system
- All accompanying system documentation is present, complete and up to date

--

Protective earth impedance (EN/IEC 62353) (Class I)	Max. reading	N/A	OK	Not OK
ME equipment / system with power connector (e.g. cold condition connector) and detachable power cord	$\leq 0.3 \Omega$ _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
ME equipment / system with non-detachable power cord	$\leq 0.3 \Omega$ _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Permanently installed ME equipment/system (e.g. ceiling mount)	$\leq 0.3 \Omega$ _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
ME system with multiple connector (e.g. table & equipment)	$\leq 0.5 \Omega$ _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Detachable power cords only	$\leq 0.1 \Omega$ _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

<b>Insulation resistance (optional)</b> (EN/IEC 62353)		<b>Reading</b>	<b>N/A</b>	<b>OK</b>	<b>Not OK</b>
Measurement at 500 VDC	$\geq 2 \text{ M}\Omega$	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Leakage current</b> (EN/IEC 62353)		<b>Max. reading</b>	<b>N/A</b>	<b>OK</b>	<b>Not OK</b>
<b>Equipment leakage current</b>					
Differential measurement:	Class I max. <b>500 <math>\mu\text{A}</math></b> Class II max. <b>100 <math>\mu\text{A}</math></b>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Direct measurement:	Class I max. <b>500 <math>\mu\text{A}</math></b> Class II max. <b>100 <math>\mu\text{A}</math></b>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Equivalent leakage current:	Class I max. <b>1000 <math>\mu\text{A}</math></b> Class II max. <b>500 <math>\mu\text{A}</math></b>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Leakage current of applied parts (BF type only)</b>					
Direct measurement:	max. <b>5000 <math>\mu\text{A}</math></b>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Function test</b> (EN/IEC 62353)			<b>N/A</b>	<b>OK</b>	<b>Not OK</b>
Function test according to SM-30-8005-116 section 9.5 incl. safety and warning functions			--	<input type="checkbox"/>	<input type="checkbox"/>
Further tested functions/parameters:			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

## Chapter 4

Test equipment used	Type	Serial number	Calibration valid until
Safety tester as per EN 62353 or VDE 0751 or EN 60601			
<b>Accessories used</b>	<b>Type</b>	<b>Serial number</b>	<b>Calibration valid until</b>

Further annexes / reports relating to the test?

Yes

No

Which:

**Additional tested equipment**

**Type**

**Serial number**

**Result**

Accessories:

_____	_____	_____
_____	_____	_____
_____	_____	_____

**Rating:**

- No safety or functional deficiencies were detected; functional tests successfully completed
- No direct risk; the detected deficiencies can be quickly corrected (see service report)
- Substantial safety deficiencies. ME equipment/system must not be operated!

Next repeat check is due in  6  12  24  36 months or not later than 08/11/2023  
(mark as applicable)

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08/11/2023

.....  
Data

.....  
Firma Service Engineer ZEISS ( Agostino Monteverde )

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# System Check Out Protocol

## OPMI Sensera, PROergo, Movena, Movena Dental, S7

APPL.	N/A		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<b>Installation</b>	→ Fill out chapter 1, 3, 4 and use either the Mobile X Customer perfect order application or the "Installation and Acceptance Report"
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<b>Preventive Maintenance</b>	→ Fill out chapter 2, 3 and 4
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<b>Corrective Maintenance</b>	→ Fill out chapter 2, 3 and 4
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<b>Safety Check</b>	→ Fill out chapter 3 and 4
		or	or
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<b>Safety Check with local report</b>	→ Fill out local Check List
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<b>Further documents attached</b>	

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	CZ Service	Customer*
<b>The Safety test according to EN/IEC 62353 is performed by:</b>	<input type="checkbox"/>	<input checked="" type="checkbox"/>

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-

Protective earth impedance (EN/IEC 62353) (Class I)	Max. reading	N/A	OK	Not OK
ME equipment / system with power connector (e.g. cold condition connector) and detachable power cord	$\leq 0.3 \Omega$ _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
ME equipment / system with non-detachable power cord	$\leq 0.3 \Omega$ _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Permanently installed ME equipment/system (e.g. ceiling mount)	$\leq 0.3 \Omega$ _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
ME system with multiple connector (e.g. table & equipment)	$\leq 0.5 \Omega$ _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Detachables power cords only	$\leq 0.1 \Omega$ _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

<b>Insulation resistance (optional)</b> (EN/IEC 62353)		<b>Reading</b>	<b>N/A</b>	<b>OK</b>	<b>Not OK</b>
Measurement at 500 VDC	<b>≥ 2 MΩ</b>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Leakage current</b> (EN/IEC 62353)		<b>Max. reading</b>	<b>N/A</b>	<b>OK</b>	<b>Not OK</b>
<b>Equipment leakage current</b>					
Differential measurement:	Class I max. <b>500 μA</b> Class II max. <b>100 μA</b>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Direct measurement:	Class I max. <b>500 μA</b> Class II max. <b>100 μA</b>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Equivalent leakage current:	Class I max. <b>1000 μA</b> Class II max. <b>500 μA</b>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Leakage current of applied parts (BF type only)</b>					
Direct measurement:	max. <b>5000 μA</b>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Function test</b> (EN/IEC 62353)			<b>N/A</b>	<b>OK</b>	<b>Not OK</b>
Function test according to SM-30-8005-116 section 9.5 incl. safety and warning functions			--	<input type="checkbox"/>	<input type="checkbox"/>
Further tested functions/parameters:			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

## Chapter 4

Test equipment used	Type	Serial number	Calibration valid until
Safety tester as per EN 62353 or VDE 0751 or EN 60601			

Accessories used	Type	Serial number	Calibration valid until

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No

Which:

Additional tested equipment	Type	Serial number	Result
Accessories:	_____	_____	_____
	_____	_____	_____
	_____	_____	_____

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